

Spetec

# Laminar Flow Systems

**Datasheet for the  
External control unit interface  
V 1.5**



# Table of Contents

|          |                               |          |   |
|----------|-------------------------------|----------|---|
| <b>1</b> | <b>INTRODUCTION</b>           | <b>3</b> |   |
| <b>2</b> | <b>CONTROL</b>                | <b>3</b> |   |
| 2.1      | Functions                     |          | 3 |
| 2.2      | Connector                     |          | 3 |
| 2.3      | Interface                     |          | 4 |
| 2.4      | Connection of In- and Outputs |          | 4 |
| <b>3</b> | <b>SPECIFICATION</b>          | <b>5</b> |   |
| 3.1      | Electrical Data               |          | 5 |
| 3.2      | Timing                        |          | 5 |
| <b>4</b> | <b>PROGRAMMING</b>            | <b>6</b> |   |
| 4.1      | Programming examples          |          | 7 |
| <b>5</b> | <b>VERSION OVERVIEW</b>       | <b>8</b> |   |

# 1 Introduction

This datasheet describes how the Spetec laminar flow systems are controlled via an external control unit (e.g. a PLC). This allows for things such as turning the system on automatically on a schedule at the beginning of production, or automatically running a flow profile without input from the operator.

The interface features potential isolation so that different ground states between the system and the external control unit have no effect.

## 2 Control

### 2.1 Functions


The following functions of the laminar flow systems can be controlled via the interface

- Power on/off
- Lights on/off
- Increase/decrease flow by one level

The following states can be read

- Power is on/off
- Lights are on/off
- Filter replacement is necessary
- A fault has occurred (fans are not turning)

### 2.2 Connector



Pinout diagram of the plug (from the solder side)

The connector is implemented as a round plug, Lumberg type WSV120.

## 2.3 Interface


Details of the interface connections

| Pin | Signal       | In/Out | Level        | Description   | Cable colour |
|-----|--------------|--------|--------------|---|--------------|
| A   | ON/OFF       | In     | high impulse | Turns the system on or off through a change in the state        | purple       |
| B   | LIGHT_ON/OFF | In     | high impulse | Turns the system lights on or off through a change in the state | green        |
| C   | FLOW+        | In     | high impulse | Increments the flow level by 1                                  | yellow       |
| D   | FLOW-        | In     | high impulse | Decrements the flow level by 1                                  | black        |
| E   | ERROR        | Out    | high         | Fans are not turning <sup>(1)</sup>                             | orange       |
|     |              |        | low          | Fans OK   |              |
| F   | FILTER       | Out    | high         | Filter replacement is necessary <sup>(2)</sup>                  | brown        |
|     |              |        | low          | Filter OK   |              |
| G   | LIGHT        | Out    | high         | The lights are turned on.                                       | white        |
|     |              |        | low          | The lights are turned off                                       |              |
| H   | POWER        | Out    | high         | The system is turned on.  | gray         |
|     |              |        | low          | The system is turned off  |              |
| J   | VCC          |        | 12-24V       | Operating voltage, 12...24V                                     | red          |
| K   | GND          |        | Ground       | Interface ground potential                                      | blue         |

### Please note:

- (1) ERROR is activated for 3-5 seconds (high) when the system is turned on. This doesn't indicate a fault as long as ERROR is deactivated after this period.
- (2) Intermittent activation indicates that the main filter is getting clogged. Multiple or permanent activation indicates that the main filter must be replaced.

## 2.4 Connection of In- and Outputs



### Please note:

- Common Ground Reference is GND (Pin K)
- Maximal voltages and levels must not exceed the electrical specification

### 3 Specification


#### 3.1 Electrical Data

|                               | Conditions            | min          | type     | max          | Unit |
|-------------------------------|-----------------------|--------------|----------|--------------|------|
| <b>General</b>                |                       |              |          |              |      |
| Operating voltage<br>$V_{CC}$ |                       | 11           | 12 or 24 | 25           | VDC  |
| Total current draw<br>(1)     | $V_{CC}=12V$          |              |          | 14           | mA   |
|                               | $V_{CC}=24V$          |              |          | 24           |      |
| <b>Inputs</b>                 |                       |              |          |              |      |
| $I_{in, high}$                |                       | 8            |          | 20           | mA   |
| $I_{in, low}$                 |                       |              |          | 0.5          | mA   |
| $U_{in, high}$                |                       | 10           |          | 25           | V    |
| $U_{in, low}$                 |                       |              |          | 1            | V    |
| <b>Outputs</b>                |                       |              |          |              |      |
| $I_{out, max}$                |                       |              |          | -20          | mA   |
| $U_{out, high}$               | $I_{out}=-1mA$        | $V_{CC}-2.3$ |          | $V_{CC}-2.0$ | V    |
|                               | $I_{out}=-10mA$       | $V_{CC}-5.4$ |          | $V_{CC}-4.9$ | V    |
| $U_{out, low}$                | $I_{out, max}=-100uA$ |              | 0,1      | 0,2          | V    |

All voltages levels are relative to the connector GND

1) Only one input activated at a time, approx. 10ma per input

#### 3.2 Timing



|                                |           | min | type | max | Unit    |
|--------------------------------|-----------|-----|------|-----|---------|
| Rise time                      | $t_r$     |     |      | 500 | $\mu s$ |
| Duration <sup>(1)</sup>        | $t_d$     | 20  |      | 100 | ms      |
| Fall time                      | $t_f$     |     |      | 500 | $\mu s$ |
| Repetition Time <sup>(2)</sup> | $t_{rep}$ | 20  |      |     | ms      |
|                                |           |     |      |     |         |

(1) The maximum duration  $t_d$  time must not be exceeded!

(2) The repetition time  $t_{rep}$  also applies to sequential control of different inputs

## 4 Programming

### During programming, please note that

- The maximum duration time  $t_d$  must not be exceeded!
- Simultaneous switching of more than one of the inputs is prohibited, since this would prevent clear detection of the system signals. All signals must be switched sequentially.
- Despite the external control unit, operation of the system via the control panel buttons and manual modification of the flow level is still possible. Therefore, a definitive flow level must be set first if the control unit is used to set a defined flow level. This is achieved by switching the FLOW+ or FLOW- input at least seven times. Then the level to be set can be selected from this state.
- When the system is turned off, all output signals are automatically set to high.

## 4.1 Programming examples

Note: CHECK(...) stands for reading the respective output

### Example 1:

Turning the system on and selecting flow level 1

```
If(CHECK(POWER) == low)
{
    ON/OFF
    7x FLOW-
}
```

### Example 2:

Changing the flow level to 5, turning the lights on

```
If(CHECK(POWER) == high)
{
    7x FLOW+
    3x FLOW-
    If(CHECK(LIGHT) = low)
    {
        LIGHT_ON/OFF
    }
}
```

### Example 3:

Turning the system on, selecting flow level 5, turning the lights off

```
If(CHECK(POWER) == low)
{
    ON/OFF
    7x FLOW+
    3x FLOW-
    If(CHECK(LIGHT) = high)
    {
        LIGHT_ON/OFF
    }
}
```

### Example 4:

Turning the lights and system off

```
If(CHECK(POWER) == high)
{
    If(CHECK(LIGHT) = high)
    {
        LIGHT_ON/OFF
    }
    ON/OFF
}
```

## 5 Version overview

| Ver. | Date       | Name, description                         |
|------|------------|---|
| 1.00 | 12.11.2012 | BachIC                                    |
| 1.10 | 09.07.2013 | Brandl, translated, 1.Freigabe            |
| 1.20 | 20.11.2014 | BachIC, some minor changes in 2.3 und 3.2 |
| 1.30 | 12.05.2015 | Kölbl, cable color expanded in 2.3        |
| 1.40 | 25.02.2016 | BachIC, connection of In/Outputs inserted |
| 1.5  | 09.07.2020 | Brandl, Programming examples adapted      |