












































PU Longlife pump tubing



















					 2 bridges	 2 bridges	 3 bridges	 3 bridges	 3 bridges
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 46.00 €	12 pieces Per package 46.00 €	12 pieces Per package 51.00 €	12 pieces Per package 51.00 €	12 pieces Per package 51.00 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number	Part number
0.381	0.0150	2.209	0.087		38-2020	381-2020	381-2220	38-2220	38-2420
0.508	0.0200	2.337	0.092		38-2025	381-2025	381-2225	38-2225	38-2425
0.635	0.0250	2.464	0.097		38-2030	381-2030	381-2230	38-2230	38-2430
0.762	0.0300	2.438	0.096		38-2032	381-2032	381-2232	38-2232	38-2432
0.889	0.0350	2.565	0.101		38-2035	381-2035	381-2235	38-2235	38-2435
1.016	0.0400	2.692	0.106		38-2040	381-2040	381-2240	38-2240	38-2440
1.143	0.0450	2.819	0.111		38-2043	381-2043	381-2243	38-2243	38-2443
1.295	0.0510	2.972	0.117		38-2046	381-2046	381-2246	38-2246	38-2446
1.422	0.0560	3.099	0.112		38-2050	381-2050	381-2250	38-2250	38-2450
1.524	0.0600	3.200	0.126		38-2052	381-2052	381-2252	38-2252	38-2452
1.651	0.0650	3.327	0.131		38-2055	381-2055	381-2255	38-2255	38-2455
1.750	0.0690	3.430	0.135		38-2058	381-2058	381-2258	38-2258	38-2458
1.854	0.0730	3.531	0.139		38-2060	381-2060	381-2260	38-2260	38-2460
2.057	0.0810	3.734	0.147		38-2065	381-2065	381-2265	38-2265	38-2465
2.286	0.0900	3.962	0.156		38-2067	381-2067	381-2267	38-2267	38-2467
2.540	0.1000	4.216	0.166		38-2070	381-2070	381-2270	38-2270	38-2470
2.794	0.1100	4.470	0.176		38-2078	381-2078	381-2278	38-2278	38-2478
3.175	0.1250	4.851	0.191		38-2080	381-2080	381-2280	38-2280	38-2480

PU Longlife – by the meter

					By the meter  3 m (118.11 inch) Per package 26.00 €	By the meter  10 m (393.70 inch) Per package 52.00 €
Inside diameter:		Outside diameter:		Colour code:		
mm	inch	mm	inch		Part number	Part number
0.381	0.0150	2.209	0.087		38-2120	38-2520
0.508	0.0200	2.337	0.092		38-2125	38-2525
0.635	0.0250	2.464	0.097		38-2130	38-2530
0.762	0.0300	2.438	0.096		38-2132	38-2532
0.889	0.0350	2.565	0.101		38-2135	38-2535
1.016	0.0400	2.692	0.106		38-2140	38-2540
1.143	0.0450	2.819	0.111		38-2143	38-2543
1.295	0.0510	2.972	0.117		38-2146	38-2546
1.422	0.0560	3.099	0.112		38-2150	38-2550
1.524	0.0600	3.200	0.126		38-2152	38-2552
1.651	0.0650	3.327	0.131		38-2155	38-2555
1.750	0.0690	3.430	0.135		38-2158	38-2558
1.854	0.0730	3.531	0.139		38-2160	38-2560
2.057	0.0810	3.734	0.147		38-2165	38-2565
2.286	0.0900	3.962	0.156		38-2167	38-2567
2.540	0.1000	4.216	0.166		38-2170	38-2570
2.794	0.1100	4.470	0.176		38-2178	38-2578
3.175	0.1250	4.851	0.191		38-2180	38-2580

PU Longlife pump tubing with "Flared End"



					 2 bridges	 2 bridges	 3 bridges	 3 bridges
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 61.00 €	12 pieces Per package 61.00 €	12 pieces Per package 66.00 €	12 pieces Per package 66.00 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number
0.381	0.0150	2.209	0.087		38-2020F	381-2020F	381-2220F	38-2220F
0.508	0.0200	2.337	0.092		38-2025F	381-2025F	381-2225F	38-2225F
0.635	0.0250	2.464	0.097		38-2030F	381-2030F	381-2230F	38-2230F
0.762	0.0300	2.438	0.096		38-2032F	381-2032F	381-2232F	38-2232F
0.889	0.0350	2.565	0.101		38-2035F	381-2035F	381-2235F	38-2235F
1.016	0.0400	2.692	0.106		38-2040F	381-2040F	381-2240F	38-2240F
1.143	0.0450	2.819	0.111		38-2043F	381-2043F	381-2243F	38-2243F
1.295	0.0510	2.972	0.117		38-2046F	381-2046F	381-2246F	38-2246F
1.422	0.0560	3.099	0.112		38-2050F	381-2050F	381-2250F	38-2250F
1.524	0.0600	3.200	0.126		38-2052F	381-2052F	381-2252F	38-2252F
1.651	0.0650	3.327	0.131		38-2055F	381-2055F	381-2255F	38-2255F