



Santoprene® TPE Pump Tubing (comparable to Pharmed/Mediprene)

Inside Diameter:		Outside Diameter:		Colour Code:	Part Number				
		mm	inch		2 bridges	2 bridges	3 bridges	3 bridges	
0.254	0.0100	2.083	0.082	orange/blue	38-6015	381-6015	381-6215	38-6215	38-6415
0.381	0.0150	2.209	0.087	orange/green	38-6020	381-6020	381-6220	38-6220	38-6420
0.508	0.0200	2.337	0.092	orange/yellow	38-6025	381-6025	381-6225	38-6225	38-6425
0.635	0.0250	2.464	0.097	orange/white	38-6030	381-6030	381-6230	38-6230	38-6430
0.762	0.0300	2.438	0.096	black/black	38-6032	381-6032	381-6232	38-6232	38-6432
0.889	0.0350	2.565	0.101	orange/orange	38-6035	381-6035	381-6235	38-6235	38-6435
1.016	0.0400	2.692	0.106	white/white	38-6040	381-6040	381-6240	38-6240	38-6440
1.143	0.0450	2.819	0.111	red/red	38-6043	381-6043	381-6243	38-6243	38-6443
1.295	0.0510	2.972	0.117	grey/grey	38-6046	381-6046	381-6246	38-6246	38-6446
1.422	0.0560	3.099	0.112	yellow/yellow	38-6050	381-6050	381-6250	38-6250	38-6450
1.524	0.0600	3.200	0.126	yellow/blue	38-6052	381-6052	381-6252	38-6252	38-6452
1.651	0.0650	3.327	0.131	blue/blue	38-6055	381-6055	381-6255	38-6255	38-6455
1.750	0.0690	3.430	0.135	blue/green	38-6058	381-6058	381-6258	38-6258	38-6458
1.854	0.0730	3.531	0.139	green/green	38-6060	381-6060	381-6260	38-6260	38-6460
2.057	0.0810	3.734	0.147	purple/purple	38-6065	381-6065	381-6265	38-6265	38-6465
2.286	0.0900	3.962	0.156	purple/black	38-6067	381-6067	381-6267	38-6267	38-6467
2.540	0.1000	4.216	0.166	purple/orange	38-6070	381-6070	381-6270	38-6270	38-6470
2.794	0.1100	4.470	0.176	purple/white	38-6078	381-6078	381-6278	38-6278	38-6478
3.175	0.1250	4.851	0.191	black/white	38-6080	381-6080	381-6280	38-6280	38-6480

Santoprene® TPE Pump Tubing (comparable to Pharmed/Mediprene) – by the meter

Inside Diameter:		Outside Diameter:		Colour Code:	Part Number	by the meter  3 m (118.11 inch) per package 32.00 €	Part Number	by the meter  10 m (393.70 inch) per package 64.00 €
mm	inch	mm	inch					
0.254	0.0100	2.083	0.082	orange/blue	38-6115		38-6515	
0.381	0.0150	2.209	0.087	orange/green	38-6120		38-6520	
0.508	0.0200	2.337	0.092	orange/yellow	38-6125		38-6525	
0.635	0.0250	2.464	0.097	orange/white	38-6130		38-6530	
0.762	0.0300	2.438	0.096	black/black	38-6132		38-6532	
0.889	0.0350	2.565	0.101	orange/orange	38-6135		38-6535	
1.016	0.0400	2.692	0.106	white/white	38-6140		38-6540	
1.143	0.0450	2.819	0.111	red/red	38-6143		38-6543	
1.295	0.0510	2.972	0.117	grey/grey	38-6146		38-6546	
1.422	0.0560	3.099	0.112	yellow/yellow	38-6150		38-6550	
1.524	0.0600	3.200	0.126	yellow/blue	38-6152		38-6552	
1.651	0.0650	3.327	0.131	blue/blue	38-6155		38-6555	
1.750	0.0690	3.430	0.135	blue/green	38-6158		38-6558	
1.854	0.0730	3.531	0.139	green/green	38-6160		38-6560	
2.057	0.0810	3.734	0.147	purple/purple	38-6165		38-6565	
2.286	0.0900	3.962	0.156	purple/black	38-6167		38-6567	
2.540	0.1000	4.216	0.166	purple/orange	38-6170		38-6570	
2.794	0.1100	4.470	0.176	purple/white	38-6178		38-6578	
3.175	0.1250	4.851	0.191	black/white	38-6180		38-6580	