

























PVC Standard – pump tubing






















Inside diameter:		Outside diameter:		Colour code:	2 bridges	2 bridges	3 bridges	3 bridges	3 bridges
mm	inch	mm	inch		12 pieces Per package 21.50 €	12 pieces Per package 21.50 €	12 pieces Per package 24.00 €	12 pieces Per package 24.00 €	12 pieces Per package 24.00 €
					Part number	Part number	Part number	Part number	Part number
0.127	0.0050	2.007	0.079		38-0005	381-0005	381-0205	38-0205	38-0405
0.191	0.0075	2.032	0.080		38-0010	381-0010	381-0210	38-0210	38-0410
0.254	0.0100	2.083	0.082		38-0015	381-0015	381-0215	38-0215	38-0415
0.381	0.0150	2.209	0.087		38-0020	381-0020	381-0220	38-0220	38-0420
0.508	0.0200	2.337	0.092		38-0025	381-0025	381-0225	38-0225	38-0425
0.635	0.0250	2.464	0.097		38-0030	381-0030	381-0230	38-0230	38-0430
0.762	0.0300	2.438	0.096		38-0032	381-0032	381-0232	38-0232	38-0432
0.889	0.0350	2.565	0.101		38-0035	381-0035	381-0235	38-0235	38-0435
1.016	0.0400	2.692	0.106		38-0040	381-0040	381-0240	38-0240	38-0440
1.143	0.0450	2.819	0.111		38-0043	381-0043	381-0243	38-0243	38-0443
1.295	0.0510	2.972	0.117		38-0046	381-0046	381-0246	38-0246	38-0446
1.422	0.0560	3.099	0.112		38-0050	381-0050	381-0250	38-0250	38-0450
1.524	0.0600	3.200	0.126		38-0052	381-0052	381-0252	38-0252	38-0452
1.651	0.0650	3.327	0.131		38-0055	381-0055	381-0255	38-0255	38-0455
1.750	0.0690	3.430	0.135		38-0058	381-0058	381-0258	38-0258	38-0458
1.854	0.0730	3.531	0.139		38-0060	381-0060	381-0260	38-0260	38-0460
2.057	0.0810	3.734	0.147		38-0065	381-0065	381-0265	38-0265	38-0465
2.286	0.0900	3.962	0.156		38-0067	381-0067	381-0267	38-0267	38-0467
2.540	0.1000	4.216	0.166		38-0070	381-0070	381-0270	38-0270	38-0470
2.794	0.1100	4.470	0.176		38-0078	381-0078	381-0278	38-0278	38-0478
3.175	0.1250	4.851	0.191		38-0080	381-0080	381-0280	38-0280	38-0480

PVC Standard – by the meter

					By the meter  3 m (118.11 inch) Per package 13.00 €	By the meter  10 m (393.70 inch) Per package 26.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number
mm	inch	mm	inch			
0.127	0.0050	2.007	0.079		38-0105	38-0505
0.191	0.0075	2.032	0.080		38-0110	38-0510
0.254	0.0100	2.083	0.082		38-0115	38-0515
0.381	0.0150	2.209	0.087		38-0120	38-0520
0.508	0.0200	2.337	0.092		38-0125	38-0525
0.635	0.0250	2.464	0.097		38-0130	38-0530
0.762	0.0300	2.438	0.096		38-0132	38-0532
0.889	0.0350	2.565	0.101		38-0135	38-0535
1.016	0.0400	2.692	0.106		38-0140	38-0540
1.143	0.0450	2.819	0.111		38-0143	38-0543
1.295	0.0510	2.972	0.117		38-0146	38-0546
1.422	0.0560	3.099	0.112		38-0150	38-0550
1.524	0.0600	3.200	0.126		38-0152	38-0552
1.651	0.0650	3.327	0.131		38-0155	38-0555
1.750	0.0690	3.430	0.135		38-0158	38-0558
1.854	0.0730	3.531	0.139		38-0160	38-0560
2.057	0.0810	3.734	0.147		38-0165	38-0565
2.286	0.0900	3.962	0.156		38-0167	38-0567
2.540	0.1000	4.216	0.166		38-0170	38-0570
2.794	0.1100	4.470	0.176		38-0178	38-0578
3.175	0.1250	4.851	0.191		38-0180	38-0580

PVC Standard pump tubing with "Flared End"



					 2 bridges	 2 bridges	 3 bridges	 3 bridges
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 36.00 €	12 pieces Per package 36.00 €	12 pieces Per package 38.50 €	12 pieces Per package 38.50 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number
0.127	0.0050	2.007	0.079		38-0005F	381-0005F	381-0205F	38-0205F
0.191	0.0075	2.032	0.080		38-0010F	381-0010F	381-0210F	38-0210F
0.254	0.0100	2.083	0.082		38-0015F	381-0015F	381-0215F	38-0215F
0.381	0.0150	2.209	0.087		38-0020F	381-0020F	381-0220F	38-0220F
0.508	0.0200	2.337	0.092		38-0025F	381-0025F	381-0225F	38-0225F
0.635	0.0250	2.464	0.097		38-0030F	381-0030F	381-0230F	38-0230F
0.762	0.0300	2.438	0.096		38-0032F	381-0032F	381-0232F	38-0232F
0.889	0.0350	2.565	0.101		38-0035F	381-0035F	381-0235F	38-0235F
1.016	0.0400	2.692	0.106		38-0040F	381-0040F	381-0240F	38-0240F
1.143	0.0450	2.819	0.111		38-0043F	381-0043F	381-0243F	38-0243F
1.295	0.0510	2.972	0.117		38-0046F	381-0046F	381-0246F	38-0246F
1.422	0.0560	3.099	0.112		38-0050F	381-0050F	381-0250F	38-0250F
1.524	0.0600	3.200	0.126		38-0052F	381-0052F	381-0252F	38-0252F
1.651	0.0650	3.327	0.131		38-0055F	381-0055F	381-0255F	38-0255F