



























# Santoprene® TPE pump tubing (comparable to Pharmed/Mediprene)



					 <b>2 bridges</b>	 <b>2 bridges</b>	 <b>3 bridges</b>	 <b>3 bridges</b>	 <b>3 bridges</b>
<b>Inside diameter:</b>		<b>Outside diameter:</b>		<b>Colour code:</b>	<b>12 pieces</b> Per package <b>56.00 €</b>	<b>12 pieces</b> Per package <b>56.00 €</b>	<b>12 pieces</b> Per package <b>60.00 €</b>	<b>12 pieces</b> Per package <b>60.00 €</b>	<b>12 pieces</b> Per package <b>60.00 €</b>
mm	inch	mm	inch		Part number	Part number	Part number	Part number	Part number
0.254	0.0100	2.083	0.082		38-6015	381-6015	381-6215	38-6215	38-6415
0.381	0.0150	2.209	0.087		38-6020	381-6020	381-6220	38-6220	38-6420
0.508	0.0200	2.337	0.092		38-6025	381-6025	381-6225	38-6225	38-6425
0.635	0.0250	2.464	0.097		38-6030	381-6030	381-6230	38-6230	38-6430
0.762	0.0300	2.438	0.096		38-6032	381-6032	381-6232	38-6232	38-6432
0.889	0.0350	2.565	0.101		38-6035	381-6035	381-6235	38-6235	38-6435
1.016	0.0400	2.692	0.106		38-6040	381-6040	381-6240	38-6240	38-6440
1.143	0.0450	2.819	0.111		38-6043	381-6043	381-6243	38-6243	38-6443
1.295	0.0510	2.972	0.117		38-6046	381-6046	381-6246	38-6246	38-6446
1.422	0.0560	3.099	0.112		38-6050	381-6050	381-6250	38-6250	38-6450
1.524	0.0600	3.200	0.126		38-6052	381-6052	381-6252	38-6252	38-6452
1.651	0.0650	3.327	0.131		38-6055	381-6055	381-6255	38-6255	38-6455
1.750	0.0690	3.430	0.135		38-6058	381-6058	381-6258	38-6258	38-6458
1.854	0.0730	3.531	0.139		38-6060	381-6060	381-6260	38-6260	38-6460
2.057	0.0810	3.734	0.147		38-6065	381-6065	381-6265	38-6265	38-6465
2.286	0.0900	3.962	0.156		38-6067	381-6067	381-6267	38-6267	38-6467
2.540	0.1000	4.216	0.166		38-6070	381-6070	381-6270	38-6270	38-6470
2.794	0.1100	4.470	0.176		38-6078	381-6078	381-6278	38-6278	38-6478
3.175	0.1250	4.851	0.191		38-6080	381-6080	381-6280	38-6280	38-6480

# Santoprene® TPE (comparable to Pharmed/Mediprene) – by the meter

					By the meter  3 m (118.11 inch) Per package 32.00 €	By the meter  10 m (393.70 inch) Per package 64.00 €
Inside diameter:		Outside diameter:		Colour code:		
mm	inch	mm	inch		Part number	Part number
0.254	0.0100	2.083	0.082		38-6115	38-6515
0.381	0.0150	2.209	0.087		38-6120	38-6520
0.508	0.0200	2.337	0.092		38-6125	38-6525
0.635	0.0250	2.464	0.097		38-6130	38-6530
0.762	0.0300	2.438	0.096		38-6132	38-6532
0.889	0.0350	2.565	0.101		38-6135	38-6535
1.016	0.0400	2.692	0.106		38-6140	38-6540
1.143	0.0450	2.819	0.111		38-6143	38-6543
1.295	0.0510	2.972	0.117		38-6146	38-6546
1.422	0.0560	3.099	0.112		38-6150	38-6550
1.524	0.0600	3.200	0.126		38-6152	38-6552
1.651	0.0650	3.327	0.131		38-6155	38-6555
1.750	0.0690	3.430	0.135		38-6158	38-6558
1.854	0.0730	3.531	0.139		38-6160	38-6560
2.057	0.0810	3.734	0.147		38-6165	38-6565
2.286	0.0900	3.962	0.156		38-6167	38-6567
2.540	0.1000	4.216	0.166		38-6170	38-6570
2.794	0.1100	4.470	0.176		38-6178	38-6578
3.175	0.1250	4.851	0.191		38-6180	38-6580