












































Solvent Flex pump tubing


















					2 bridges  12 pieces Per package 37.00 €	2 bridges  12 pieces Per package 37.00 €	3 bridges  12 pieces Per package 41.00 €	3 bridges  12 pieces Per package 41.00 €	3 bridges  12 pieces Per package 41.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number	Part number	Part number	Part number
mm	inch	mm	inch						
0.381	0.0150	2.209	0.087		38-4020	381-4020	381-4220	38-4220	38-4420
0.508	0.0200	2.337	0.092		38-4025	381-4025	381-4225	38-4225	38-4425
0.635	0.0250	2.464	0.097		38-4030	381-4030	381-4230	38-4230	38-4430
0.762	0.0300	2.438	0.096		38-4032	381-4032	381-4232	38-4232	38-4432
0.889	0.0350	2.565	0.101		38-4035	381-4035	381-4235	38-4235	38-4435
1.016	0.0400	2.692	0.106		38-4040	381-4040	381-4240	38-4240	38-4440
1.143	0.0450	2.819	0.111		38-4043	381-4043	381-4243	38-4243	38-4443
1.295	0.0510	2.972	0.117		38-4046	381-4046	381-4246	38-4246	38-4446
1.422	0.0560	3.099	0.112		38-4050	381-4050	381-4250	38-4250	38-4450
1.524	0.0600	3.200	0.126		38-4052	381-4052	381-4252	38-4252	38-4452
1.651	0.0650	3.327	0.131		38-4055	381-4055	381-4255	38-4255	38-4455
1.750	0.0690	3.430	0.135		38-4058	381-4058	381-4258	38-4258	38-4458
1.854	0.0730	3.531	0.139		38-4060	381-4060	381-4260	38-4260	38-4460
2.057	0.0810	3.734	0.147		38-4065	381-4065	381-4265	38-4265	38-4465
2.286	0.0900	3.962	0.156		38-4067	381-4067	381-4267	38-4267	38-4467
2.540	0.1000	4.216	0.166		38-4070	381-4070	381-4270	38-4270	38-4470
2.794	0.1100	4.470	0.176		38-4078	381-4078	381-4278	38-4278	38-4478
3.175	0.1250	4.851	0.191		38-4080	381-4080	381-4280	38-4280	38-4480

Solvent Flex – by the meter

					By the meter  3 m (118.11 inch) Per package 21.00 €	By the meter  10 m (393.70 inch) Per package 42.00 €
Inside diameter:		Outside diameter:		Colour code:		
mm	inch	mm	inch		Part number	Part number
0.381	0.0150	2.209	0.087		38-4120	38-4520
0.508	0.0200	2.337	0.092		38-4125	38-4525
0.635	0.0250	2.464	0.097		38-4130	38-4530
0.762	0.0300	2.438	0.096		38-4132	38-4532
0.889	0.0350	2.565	0.101		38-4135	38-4535
1.016	0.0400	2.692	0.106		38-4140	38-4540
1.143	0.0450	2.819	0.111		38-4143	38-4543
1.295	0.0510	2.972	0.117		38-4146	38-4546
1.422	0.0560	3.099	0.112		38-4150	38-4550
1.524	0.0600	3.200	0.126		38-4152	38-4552
1.651	0.0650	3.327	0.131		38-4155	38-4555
1.750	0.0690	3.430	0.135		38-4158	38-4558
1.854	0.0730	3.531	0.139		38-4160	38-4560
2.057	0.0810	3.734	0.147		38-4165	38-4565
2.286	0.0900	3.962	0.156		38-4167	38-4567
2.540	0.1000	4.216	0.166		38-4170	38-4570
2.794	0.1100	4.470	0.176		38-4178	38-4578
3.175	0.1250	4.851	0.191		38-4180	38-4580

Solvent Flex pump tubing with "Flared End"

					 2 bridges	 2 bridges	 3 bridges	 3 bridges
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 52.00 €	12 pieces Per package 52.00 €	12 pieces Per package 56.00 €	12 pieces Per package 56.00 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number
0.381	0.0150	2.209	0.087		38-4020F	381-4020F	381-4220F	38-4220F
0.508	0.0200	2.337	0.092		38-4025F	381-4025F	381-4225F	38-4225F
0.635	0.0250	2.464	0.097		38-4030F	381-4030F	381-4230F	38-4230F
0.762	0.0300	2.438	0.096		38-4032F	381-4032F	381-4232F	38-4232F
0.889	0.0350	2.565	0.101		38-4035F	381-4035F	381-4235F	38-4235F
1.016	0.0400	2.692	0.106		38-4040F	381-4040F	381-4240F	38-4240F
1.143	0.0450	2.819	0.111		38-4043F	381-4043F	381-4243F	38-4243F
1.295	0.0510	2.972	0.117		38-4046F	381-4046F	381-4246F	38-4246F
1.422	0.0560	3.099	0.112		38-4050F	381-4050F	381-4250F	38-4250F
1.524	0.0600	3.200	0.126		38-4052F	381-4052F	381-4252F	38-4252F
1.651	0.0650	3.327	0.131		38-4055F	381-4055F	381-4255F	38-4255F