



Tubing

Directly from the manufacturer



Variety is where our strength



Spetec tubing is available a number of different materials: PVC, polyurethane, fluorinated rubber comparable to Viton[®], Santoprene[®] and silicone. The tubing is available in different configurations for different applications.

Our entire range covers more than 450 items and includes tubing with an internal diameter of 0.12 through 3.2 mm, which we supply to end customers, retailers and equipment manufacturers.



Do you need a solution tailored to your requirements? Just send us the specifications you require and we will work alongside you to find an optimum solution with short delivery times!

And we are happy to provide you with custom tubing. All quality grades and configurations are available at short notice.



Tubing direct from the manufacturer

	Page		Page
The materials	4	Fluorinated rubber –	
The sizes	6	Pump tubing	17
Tubing		By the meter	18
PVC Standard –		Dual pump tubing	19
Pump tubing	7	Santoprene® –	
By the meter	8	Pump tubing	20
Flared End	9	By the meter	21
Dual pump tubing	10	Dual pump tubing	22
Solvent Flex –		Silicone –	
Pump tubing	11	Pump tubing	23
By the meter	12	By the meter	24
Flared End	13	Dual pump tubing	25
PU Longlife –		Tubing connectors	26
Pump tubing	14	Special tubing	27
By the meter	15	Single-element and multi-element standards	28
Dual pump tubing	16	OEM and standalone pumps	29

All prices are subject to VAT.
Errors and omissions excepted. Prices subject to change without notice.

Tubing – The materials



Please also request our data sheet!

PVC Standard

Thermoplastic polymer that is transformed into a soft, multipurpose tubing material by the addition of plasticizers and stabilizers.

Service life	300 operating hours
Material	Polyvinyl chloride
Specific gravity	1.20 g/cm ³
Hardness	68 ± 5 Shore A
Temperature range for pump operation	18 – 28 °C
Temperature range of raw material	5 – 40 °C
Color	Transparent
Overprinting	Internal diameter in inches – S
Chemical resistance	Dilute alkalis, dilute acids, ammonia
Comparable to	–

Solvent Flex

Thermoplastic polymer offering enhanced strength thanks to the use of special additives.

Service life	200 operating hours
Material	Polyvinyl chloride
Specific gravity	1.18g/cm ³
Hardness	65 ± 5 Shore A
Temperature range for pump operation	18 – 28 °C
Temperature range of raw material	5 – 40 °C
Color	Yellow – transparent
Overprinting	None
Chemical resistance	Oils
Comparable to	–

PU Longlife

Thermoplastic polyether polyurethane with excellent resistance to hydrolysis.

Service life	500 operating hours
Material	Polyurethane
Specific gravity	1.14 g/cm ³
Hardness	75 ± 5 Shore A
Temperature range for pump operation	18 – 28 °C
Temperature range of raw material	5 – 40 °C
Color	Transparent
Overprinting	None
Chemical resistance	Alcohols, benzene, diesel and kerosene
Comparable to	–

Fluorinated rubber comparable to Viton®

In its elastic form, fluorinated rubber is characterized by its excellent chemical and thermal properties.

Service life	150 operating hours
Material	Fluoroelastomer
Specific gravity	1.89 g/cm ³
Hardness	60 ± 5 Shore A
Temperature range for pump operation	18 – 28 °C
Temperature range of raw material	-10 – 200°C
Color	Black
Overprinting	None
Chemical resistance	Acids, alkalis, oils
Comparable to	Viton®

Santoprene® comparable to Pharmed®/ Mediprene®

As a thermoplastic elastomer, Santoprene® provides excellent chemical properties and mechanical strength.

Service life	300 operating hours
Material	Thermoplastic vulcanized material
Specific gravity	0.97g/cm ³
Hardness	64 ± 5 Shore A
Temperature range for pump operation	18 – 28 °C
Temperature range of raw material	-40 – 125°C
Color	Beige
Overprinting	None
Chemical resistance	Dilute acids, solvents, acetone
Comparable to	Mediprene®, Pharmed®

Silicone


As an autoclave-safe material, this silicone-based tubing is particularly suitable for use in the foodstuffs industry.

Service life	300 operating hours
Material	Silicone
Specific gravity	1.14 g/cm ³
Hardness	60 ± 5 Shore A
Temperature range for pump operation	18 – 28 °C
Temperature range of raw material	-40 – 200°C
Color	Translucent
Overprinting	None
Chemical resistance	Dilute acids
Comparable to	–


Tubing – The sizes




2 bridges




2 bridges




3 bridges




3 bridges




By the meter: 3 m or 10 m









2 bridges



PVC Standard pump tubing







					 2 bridges	 2 bridges	 3 bridges	 3 bridges	 3 bridges
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 22.50 €	12 pieces Per package 22.50 €	12 pieces Per package 25.00 €	12 pieces Per package 25.00 €	12 pieces Per package 25.00 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number	Part number
0,191	0,0075	2,032	0,080		38-0010	381-0010	381-0210	38-0210	38-0410
0,254	0,0100	2,083	0,082		38-0015	381-0015	381-0215	38-0215	38-0415
0,381	0,0150	2,209	0,087		38-0020	381-0020	381-0220	38-0220	38-0420
0,440	0,0170	2,261	0,089		38-0022	381-0022	381-0222	38-0222	38-0422
0,508	0,0200	2,337	0,092		38-0025	381-0025	381-0225	38-0225	38-0425
0,635	0,0250	2,438	0,097		38-0030	381-0030	381-0230	38-0230	38-0430
0,762	0,0300	2,464	0,096		38-0032	381-0032	381-0232	38-0232	38-0432
0,889	0,0350	2,565	0,101		38-0035	381-0035	381-0235	38-0235	38-0435
1,016	0,0400	2,692	0,106		38-0040	381-0040	381-0240	38-0240	38-0440
1,143	0,0450	2,819	0,111		38-0043	381-0043	381-0243	38-0243	38-0443
1,295	0,0510	2,972	0,117		38-0046	381-0046	381-0246	38-0246	38-0446
1,422	0,0560	3,099	0,112		38-0050	381-0050	381-0250	38-0250	38-0450
1,524	0,0600	3,200	0,126		38-0052	381-0052	381-0252	38-0252	38-0452
1,651	0,0650	3,327	0,131		38-0055	381-0055	381-0255	38-0255	38-0455
1,750	0,0690	3,430	0,135		38-0058	381-0058	381-0258	38-0258	38-0458
1,854	0,0730	3,531	0,139		38-0060	381-0060	381-0260	38-0260	38-0460
2,057	0,0810	3,734	0,147		38-0065	381-0065	381-0265	38-0265	38-0465
2,286	0,0900	3,962	0,156		38-0067	381-0067	381-0267	38-0267	38-0467
2,540	0,1000	4,216	0,166		38-0070	381-0070	381-0270	38-0270	38-0470
2,794	0,1100	4,470	0,176		38-0078	381-0078	381-0278	38-0278	38-0478
3,175	0,1250	4,851	0,191		38-0080	381-0080	381-0280	38-0280	38-0480

PVC Standard – by the meter


					By the meter  3 m (118.11 inch) Per package 13.50 €	By the meter  10 m (393.70 inch) Per package 27.00€
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number
mm	inch	mm	inch			
0,191	0,0075	2,032	0,080		38-0110	38-0510
0,254	0,0100	2,083	0,082		38-0115	38-0515
0,381	0,0150	2,209	0,087		38-0120	38-0520
0,440	0,0170	2,261	0,089		38-0122	38-0522
0,508	0,0200	2,337	0,092		38-0125	38-0525
0,635	0,0250	2,438	0,097		38-0130	38-0530
0,762	0,0300	2,464	0,096		38-0132	38-0532
0,889	0,0350	2,565	0,101		38-0135	38-0535
1,016	0,0400	2,692	0,106		38-0140	38-0540
1,143	0,0450	2,819	0,111		38-0143	38-0543
1,295	0,0510	2,972	0,117		38-0146	38-0546
1,422	0,0560	3,099	0,112		38-0150	38-0550
1,524	0,0600	3,200	0,126		38-0152	38-0552
1,651	0,0650	3,327	0,131		38-0155	38-0555
1,750	0,0690	3,430	0,135		38-0158	38-0558
1,854	0,0730	3,531	0,139		38-0160	38-0560
2,057	0,0810	3,734	0,147		38-0165	38-0565
2,286	0,0900	3,962	0,156		38-0167	38-0567
2,540	0,1000	4,216	0,166		38-0170	38-0570
2,794	0,1100	4,470	0,176		38-0178	38-0578
3,175	0,1250	4,851	0,191		38-0180	38-0580

PVC Standard pump tubing with "Flared End"








					2 bridges  12 pieces Per package 37.50 €	2 bridges  12 pieces Per package 37.50 €	3 bridges  12 pieces Per package 40.00 €	3 bridges  12 pieces Per package 40.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number	Part number	Part number
mm	inch	mm	inch					
0,191	0,0075	2,032	0,080		38-0010F	381-0010F	381-0210F	38-0210F
0,254	0,0100	2,083	0,082		38-0015F	381-0015F	381-0215F	38-0215F
0,381	0,0150	2,209	0,087		38-0020F	381-0020F	381-0220F	38-0220F
0,440	0,0170	2,261	0,089		38-0022F	381-0022F	381-0222F	38-0222F
0,508	0,0200	2,337	0,092		38-0025F	381-0025F	381-0225F	38-0225F
0,635	0,0250	2,438	0,097		38-0030F	381-0030F	381-0230F	38-0230F
0,762	0,0300	2,464	0,096		38-0032F	381-0032F	381-0232F	38-0232F
0,889	0,0350	2,565	0,101		38-0035F	381-0035F	381-0235F	38-0235F
1,016	0,0400	2,692	0,106		38-0040F	381-0040F	381-0240F	38-0240F
1,143	0,0450	2,819	0,111		38-0043F	381-0043F	381-0243F	38-0243F
1,295	0,0510	2,972	0,117		38-0046F	381-0046F	381-0246F	38-0246F
1,422	0,0560	3,099	0,112		38-0050F	381-0050F	381-0250F	38-0250F
1,524	0,0600	3,200	0,126		38-0052F	381-0052F	381-0252F	38-0252F
1,651	0,0650	3,327	0,131		38-0055F	381-0055F	381-0255F	38-0255F

PVC Standard double pump tubing


Pairs of tubing connected with Y-pieces on both sides extended with a tube of 25cm					 2 bridges 152 mm 250 mm 406 mm 250 mm per package 6 pieces 75.00 €
Inside diameter:		Outside diameter:		Colour code:	
mm	inch	mm	inch		Part number
1.016	0.0400	2.692	0.106		38-0340
1.143	0.0450	2.819	0.111		38-0343
1.295	0.0510	2.972	0.117		38-0346
1.422	0.0560	3.099	0.112		38-0350
1.524	0.0600	3.200	0.126		38-0352
1.651	0.0650	3.327	0.131		38-0355
1.750	0.0690	3.430	0.135		38-0358
1.854	0.0730	3.531	0.139		38-0360
2.057	0.0810	3.734	0.147		38-0365
2.286	0.0900	3.962	0.156		38-0367
2.540	0.1000	4.216	0.166		38-0370
2.794	0.1100	4.470	0.176		38-0378
3.175	0.1250	4.851	0.191		38-0380

Solvent Flex pump tubing







					2 bridges 	2 bridges 	3 bridges 	3 bridges 	3 bridges 
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 38.50 €	12 pieces Per package 38.50 €	12 pieces Per package 42.50 €	12 pieces Per package 42.50 €	12 pieces Per package 42.50 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number	Part number
0.381	0.0150	2.209	0.087		38-4020	381-4020	381-4220	38-4220	38-4420
0.508	0.0200	2.337	0.092		38-4025	381-4025	381-4225	38-4225	38-4425
0.635	0.0250	2.438	0.097		38-4030	381-4030	381-4230	38-4230	38-4430
0.762	0.0300	2.464	0.096		38-4032	381-4032	381-4232	38-4232	38-4432
0.889	0.0350	2.565	0.101		38-4035	381-4035	381-4235	38-4235	38-4435
1.016	0.0400	2.692	0.106		38-4040	381-4040	381-4240	38-4240	38-4440
1.143	0.0450	2.819	0.111		38-4043	381-4043	381-4243	38-4243	38-4443
1.295	0.0510	2.972	0.117		38-4046	381-4046	381-4246	38-4246	38-4446
1.422	0.0560	3.099	0.112		38-4050	381-4050	381-4250	38-4250	38-4450
1.524	0.0600	3.200	0.126		38-4052	381-4052	381-4252	38-4252	38-4452
1.651	0.0650	3.327	0.131		38-4055	381-4055	381-4255	38-4255	38-4455
1.750	0.0690	3.430	0.135		38-4058	381-4058	381-4258	38-4258	38-4458
1.854	0.0730	3.531	0.139		38-4060	381-4060	381-4260	38-4260	38-4460
2.057	0.0810	3.734	0.147		38-4065	381-4065	381-4265	38-4265	38-4465
2.286	0.0900	3.962	0.156		38-4067	381-4067	381-4267	38-4267	38-4467
2.540	0.1000	4.216	0.166		38-4070	381-4070	381-4270	38-4270	38-4470
2.794	0.1100	4.470	0.176		38-4078	381-4078	381-4278	38-4278	38-4478
3.175	0.1250	4.851	0.191		38-4080	381-4080	381-4280	38-4280	38-4480

Solvent Flex – by the meter

					By the meter  3 m (118.11 inch) Per package 22.00 €	By the meter  10 m (393.70 inch) Per package 44.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number
mm	inch	mm	inch			
0.381	0.0150	2.209	0.087		38-4120	38-4520
0.508	0.0200	2.337	0.092		38-4125	38-4525
0.635	0.0250	2.438	0.097		38-4130	38-4530
0.762	0.0300	2.464	0.096		38-4132	38-4532
0.889	0.0350	2.565	0.101		38-4135	38-4535
1.016	0.0400	2.692	0.106		38-4140	38-4540
1.143	0.0450	2.819	0.111		38-4143	38-4543
1.295	0.0510	2.972	0.117		38-4146	38-4546
1.422	0.0560	3.099	0.112		38-4150	38-4550
1.524	0.0600	3.200	0.126		38-4152	38-4552
1.651	0.0650	3.327	0.131		38-4155	38-4555
1.750	0.0690	3.430	0.135		38-4158	38-4558
1.854	0.0730	3.531	0.139		38-4160	38-4560
2.057	0.0810	3.734	0.147		38-4165	38-4565
2.286	0.0900	3.962	0.156		38-4167	38-4567
2.540	0.1000	4.216	0.166		38-4170	38-4570
2.794	0.1100	4.470	0.176		38-4178	38-4578
3.175	0.1250	4.851	0.191		38-4180	38-4580

Solvent Flex pump tubing with “Flared End”



					2 bridges  12 pieces Per package 54.00 €	2 bridges  12 pieces Per package 54.00 €	3 bridges  12 pieces Per package 58.50 €	3 bridges  12 pieces Per package 58.50 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number	Part number	Part number
mm	inch	mm	inch					
0.381	0.0150	2.209	0.087		38-4020F	381-4020F	381-4220F	38-4220F
0.508	0.0200	2.337	0.092		38-4025F	381-4025F	381-4225F	38-4225F
0.635	0.0250	2.438	0.097		38-4030F	381-4030F	381-4230F	38-4230F
0.762	0.0300	2.464	0.096		38-4032F	381-4032F	381-4232F	38-4232F
0.889	0.0350	2.565	0.101		38-4035F	381-4035F	381-4235F	38-4235F
1.016	0.0400	2.692	0.106		38-4040F	381-4040F	381-4240F	38-4240F
1.143	0.0450	2.819	0.111		38-4043F	381-4043F	381-4243F	38-4243F
1.295	0.0510	2.972	0.117		38-4046F	381-4046F	381-4246F	38-4246F
1.422	0.0560	3.099	0.112		38-4050F	381-4050F	381-4250F	38-4250F
1.524	0.0600	3.200	0.126		38-4052F	381-4052F	381-4252F	38-4252F
1.651	0.0650	3.327	0.131		38-4055F	381-4055F	381-4255F	38-4255F

PU Longlife pump tubing




					 2 bridges 12 pieces Per package 46.00 €	 2 bridges 12 pieces Per package 46.00 €	 3 bridges 12 pieces Per package 51.00 €	 3 bridges 12 pieces Per package 51.00 €	 3 bridges 12 pieces Per package 51.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number	Part number	Part number	Part number
mm	inch	mm	inch						
0.381	0.0150	2.209	0.087		38-2020	381-2020	381-2220	38-2220	38-2420
0.508	0.0200	2.337	0.092		38-2025	381-2025	381-2225	38-2225	38-2425
0.635	0.0250	2.438	0.097		38-2030	381-2030	381-2230	38-2230	38-2430
0.762	0.0300	2.464	0.096		38-2032	381-2032	381-2232	38-2232	38-2432
0.889	0.0350	2.565	0.101		38-2035	381-2035	381-2235	38-2235	38-2435
1.016	0.0400	2.692	0.106		38-2040	381-2040	381-2240	38-2240	38-2440
1.143	0.0450	2.819	0.111		38-2043	381-2043	381-2243	38-2243	38-2443
1.295	0.0510	2.972	0.117		38-2046	381-2046	381-2246	38-2246	38-2446
1.422	0.0560	3.099	0.112		38-2050	381-2050	381-2250	38-2250	38-2450
1.524	0.0600	3.200	0.126		38-2052	381-2052	381-2252	38-2252	38-2452
1.651	0.0650	3.327	0.131		38-2055	381-2055	381-2255	38-2255	38-2455
1.750	0.0690	3.430	0.135		38-2058	381-2058	381-2258	38-2258	38-2458
1.854	0.0730	3.531	0.139		38-2060	381-2060	381-2260	38-2260	38-2460
2.057	0.0810	3.734	0.147		38-2065	381-2065	381-2265	38-2265	38-2465
2.286	0.0900	3.962	0.156		38-2067	381-2067	381-2267	38-2267	38-2467
2.540	0.1000	4.216	0.166		38-2070	381-2070	381-2270	38-2270	38-2470
2.794	0.1100	4.470	0.176		38-2078	381-2078	381-2278	38-2278	38-2478
3.175	0.1250	4.851	0.191		38-2080	381-2080	381-2280	38-2280	38-2480

PU Longlife – by the meter






					By the meter  3 m (118.11 inch) Per package 26.00 €	By the meter  10 m (393.70 inch) Per package 52.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number
mm	inch	mm	inch			
0.381	0.0150	2.209	0.087		38-2120	38-2520
0.508	0.0200	2.337	0.092		38-2125	38-2525
0.635	0.0250	2.438	0.097		38-2130	38-2530
0.762	0.0300	2.464	0.096		38-2132	38-2532
0.889	0.0350	2.565	0.101		38-2135	38-2535
1.016	0.0400	2.692	0.106		38-2140	38-2540
1.143	0.0450	2.819	0.111		38-2143	38-2543
1.295	0.0510	2.972	0.117		38-2146	38-2546
1.422	0.0560	3.099	0.112		38-2150	38-2550
1.524	0.0600	3.200	0.126		38-2152	38-2552
1.651	0.0650	3.327	0.131		38-2155	38-2555
1.750	0.0690	3.430	0.135		38-2158	38-2558
1.854	0.0730	3.531	0.139		38-2160	38-2560
2.057	0.0810	3.734	0.147		38-2165	38-2565
2.286	0.0900	3.962	0.156		38-2167	38-2567
2.540	0.1000	4.216	0.166		38-2170	38-2570
2.794	0.1100	4.470	0.176		38-2178	38-2578
3.175	0.1250	4.851	0.191		38-2180	38-2580

PU longlife double pump tubing

Pairs of tubing connected with Y-pieces on both sides extended with a tube of 25cm					 2 bridges 152 mm 250 mm 406 mm 250 mm per package 6 pieces 89.00 €
Inside diameter:		Outside diameter:		Colour code:	
mm	inch	mm	inch		Part number
1.016	0.0400	2.692	0.106		38-2340
1.143	0.0450	2.819	0.111		38-2343
1.295	0.0510	2.972	0.117		38-2346
1.422	0.0560	3.099	0.112		38-2350
1.524	0.0600	3.200	0.126		38-2352
1.651	0.0650	3.327	0.131		38-2355
1.750	0.0690	3.430	0.135		38-2358
1.854	0.0730	3.531	0.139		38-2360
2.057	0.0810	3.734	0.147		38-2365
2.286	0.0900	3.962	0.156		38-2367
2.540	0.1000	4.216	0.166		38-2370
2.794	0.1100	4.470	0.176		38-2378
3.175	0.1250	4.851	0.191		38-2380

Fluoro rubber pump tubing (comparable to Viton®)




					2 bridges 	2 bridges 	3 bridges 	3 bridges 	3 bridges 
Inside diameter:		Outside diameter:		Colour code:	12 pieces Per package 115.00 €	12 pieces Per package 115.00 €	12 pieces Per package 132.00 €	12 pieces Per package 132.00 €	12 pieces Per package 132.00 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number	Part number
0.635	0.0250	2.438	0.097		38-5630	381-5630	381-5230	38-5230	38-5430
0.762	0.0300	2.464	0.096		38-5632	381-5632	381-5232	38-5232	38-5432
0.889	0.0350	2.565	0.101		38-5635	381-5635	381-5235	38-5235	38-5435
1.016	0.0400	2.692	0.106		38-5640	381-5640	381-5240	38-5240	38-5440
1.143	0.0450	2.819	0.111		38-5643	381-5643	381-5243	38-5243	38-5443
1.295	0.0510	2.972	0.117		38-5646	381-5646	381-5246	38-5246	38-5446
1.422	0.0560	3.099	0.112		38-5650	381-5650	381-5250	38-5250	38-5450
1.524	0.0600	3.200	0.126		38-5652	381-5652	381-5252	38-5252	38-5452
1.651	0.0650	3.327	0.131		38-5655	381-5655	381-5255	38-5255	38-5455
1.750	0.0690	3.430	0.135		38-5658	381-5658	381-5258	38-5258	38-5458
1.854	0.0730	3.531	0.139		38-5660	381-5660	381-5260	38-5260	38-5460
2.057	0.0810	3.734	0.147		38-5665	381-5665	381-5265	38-5265	38-5465
2.286	0.0900	3.962	0.156		38-5667	381-5667	381-5267	38-5267	38-5467
2.540	0.1000	4.216	0.166		38-5670	381-5670	381-5270	38-5270	38-5470
2.794	0.1100	4.470	0.176		38-5678	381-5678	381-5278	38-5278	38-5478
3.175	0.1250	4.851	0.191		38-5680	381-5680	381-5280	38-5280	38-5480

Fluoro rubber – by the meter (comparable to Viton®)

					By the meter  3 m (118.11 inch) Per package 69.00 €	By the meter  10 m (393.70 inch) Per package 138.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number
mm	inch	mm	inch			
0.635	0.0250	2.438	0.097		38-5130	38-5530
0.762	0.0300	2.464	0.096		38-5132	38-5532
0.889	0.0350	2.565	0.101		38-5135	38-5535
1.016	0.0400	2.692	0.106		38-5140	38-5540
1.143	0.0450	2.819	0.111		38-5143	38-5543
1.295	0.0510	2.972	0.117		38-5146	38-5546
1.422	0.0560	3.099	0.112		38-5150	38-5550
1.524	0.0600	3.200	0.126		38-5152	38-5552
1.651	0.0650	3.327	0.131		38-5155	38-5555
1.750	0.0690	3.430	0.135		38-5158	38-5558
1.854	0.0730	3.531	0.139		38-5160	38-5560
2.057	0.0810	3.734	0.147		38-5165	38-5565
2.286	0.0900	3.962	0.156		38-5167	38-5567
2.540	0.1000	4.216	0.166		38-5170	38-5570
2.794	0.1100	4.470	0.176		38-5178	38-5578
3.175	0.1250	4.851	0.191		38-5180	38-5580






Fluoro rubber double pump tubing (comparable to Viton®)



Pairs of tubing connected with Y-pieces on both sides extended with a tube of 25cm					2 bridges 250 mm 406 mm 250 mm per package 6 pieces 217.50 €
Inside diameter:		Outside diameter:		Colour code:	Part number
mm	inch	mm	inch		
1.016	0.0400	2.692	0.106		38-5340
1.143	0.0450	2.819	0.111		38-5343
1.295	0.0510	2.972	0.117		38-5346
1.422	0.0560	3.099	0.112		38-5350
1.524	0.0600	3.200	0.126		38-5352
1.651	0.0650	3.327	0.131		38-5355
1.750	0.0690	3.430	0.135		38-5358
1.854	0.0730	3.531	0.139		38-5360
2.057	0.0810	3.734	0.147		38-5365
2.286	0.0900	3.962	0.156		38-5367
2.540	0.1000	4.216	0.166		38-5370
2.794	0.1100	4.470	0.176		38-5378
3.175	0.1250	4.851	0.191		38-5380

Santoprene® pump tubing (comparable to Pharmed®/Mediprene®)




					 2 bridges 12 pieces Per package 62.00 €	 2 bridges 12 pieces Per package 62.00 €	 3 bridges 12 pieces Per package 66.00 €	 3 bridges 12 pieces Per package 66.00 €	 3 bridges 12 pieces Per package 66.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number	Part number	Part number	Part number
mm	inch	mm	inch						
0.254	0.0100	2.083	0.082		38-6015	381-6015	381-6215	38-6215	38-6415
0.381	0.0150	2.209	0.087		38-6020	381-6020	381-6220	38-6220	38-6420
0.508	0.0200	2.337	0.092		38-6025	381-6025	381-6225	38-6225	38-6425
0.635	0.0250	2.438	0.097		38-6030	381-6030	381-6230	38-6230	38-6430
0.762	0.0300	2.464	0.096		38-6032	381-6032	381-6232	38-6232	38-6432
0.889	0.0350	2.565	0.101		38-6035	381-6035	381-6235	38-6235	38-6435
1.016	0.0400	2.692	0.106		38-6040	381-6040	381-6240	38-6240	38-6440
1.143	0.0450	2.819	0.111		38-6043	381-6043	381-6243	38-6243	38-6443
1.295	0.0510	2.972	0.117		38-6046	381-6046	381-6246	38-6246	38-6446
1.422	0.0560	3.099	0.112		38-6050	381-6050	381-6250	38-6250	38-6450
1.524	0.0600	3.200	0.126		38-6052	381-6052	381-6252	38-6252	38-6452
1.651	0.0650	3.327	0.131		38-6055	381-6055	381-6255	38-6255	38-6455
1.750	0.0690	3.430	0.135		38-6058	381-6058	381-6258	38-6258	38-6458
1.854	0.0730	3.531	0.139		38-6060	381-6060	381-6260	38-6260	38-6460
2.057	0.0810	3.734	0.147		38-6065	381-6065	381-6265	38-6265	38-6465
2.286	0.0900	3.962	0.156		38-6067	381-6067	381-6267	38-6267	38-6467
2.540	0.1000	4.216	0.166		38-6070	381-6070	381-6270	38-6270	38-6470
2.794	0.1100	4.470	0.176		38-6078	381-6078	381-6278	38-6278	38-6478
3.175	0.1250	4.851	0.191		38-6080	381-6080	381-6280	38-6280	38-6480

Santoprene® – by the meter


(comparable to Pharmed®/Mediprene®)




					By the meter  3 m (118.11 inch) Per package 33.50 €	By the meter  10 m (393.70 inch) Per package 67.00 €
Inside diameter:		Outside diameter:		Colour code:	Part number	Part number
mm	inch	mm	inch			
0.254	0.0100	2.083	0.082		38-6115	38-6515
0.381	0.0150	2.209	0.087		38-6120	38-6520
0.508	0.0200	2.337	0.092		38-6125	38-6525
0.635	0.0250	2.438	0.097		38-6130	38-6530
0.762	0.0300	2.464	0.096		38-6132	38-6532
0.889	0.0350	2.565	0.101		38-6135	38-6535
1.016	0.0400	2.692	0.106		38-6140	38-6540
1.143	0.0450	2.819	0.111		38-6143	38-6543
1.295	0.0510	2.972	0.117		38-6146	38-6546
1.422	0.0560	3.099	0.112		38-6150	38-6550
1.524	0.0600	3.200	0.126		38-6152	38-6552
1.651	0.0650	3.327	0.131		38-6155	38-6555
1.750	0.0690	3.430	0.135		38-6158	38-6558
1.854	0.0730	3.531	0.139		38-6160	38-6560
2.057	0.0810	3.734	0.147		38-6165	38-6565
2.286	0.0900	3.962	0.156		38-6167	38-6567
2.540	0.1000	4.216	0.166		38-6170	38-6570
2.794	0.1100	4.470	0.176		38-6178	38-6578
3.175	0.1250	4.851	0.191		38-6180	38-6580










Santoprene® double pump tubing

(comparable to Pharmed®/Mediprene®)

Pairs of tubing connected with Y-pieces on both sides extended with a tube of 25cm					 2 bridges 152 mm 250 mm 406 mm 250 mm per package 6 pieces 125.00 €
Inside diameter:		Outside diameter:		Colour code:	
mm	inch	mm	inch		Part number
1.016	0.0400	2.692	0.106		38-6340
1.143	0.0450	2.819	0.111		38-6343
1.295	0.0510	2.972	0.117		38-6346
1.422	0.0560	3.099	0.112		38-6350
1.524	0.0600	3.200	0.126		38-6352
1.651	0.0650	3.327	0.131		38-6355
1.750	0.0690	3.430	0.135		38-6358
1.854	0.0730	3.531	0.139		38-6360
2.057	0.0810	3.734	0.147		38-6365
2.286	0.0900	3.962	0.156		38-6367
2.540	0.1000	4.216	0.166		38-6370
2.794	0.1100	4.470	0.176		38-6378
3.175	0.1250	4.851	0.191		38-6380

Silicone pump tubing




					 2 bridges	 2 bridges	 3 bridges	 3 bridges	 3 bridges
Inside diameter:		Outside diameter:		Colour code:	6 pieces Per package 44.00 €	6 pieces Per package 44.00 €	6 pieces Per package 48.00 €	6 pieces Per package 48.00 €	6 pieces Per package 48.00 €
mm	inch	mm	inch		Part number	Part number	Part number	Part number	Part number
0.508	0.0200	2.337	0.092		38-3025	381-3025	381-3225	38-3225	38-3425
0.635	0.0250	2.438	0.097		38-3030	381-3030	381-3230	38-3230	38-3430
0.762	0.0300	2.464	0.096		38-3032	381-3032	381-3232	38-3232	38-3432
0.889	0.0350	2.565	0.101		38-3035	381-3035	381-3235	38-3235	38-3435
1.016	0.0400	2.692	0.106		38-3040	381-3040	381-3240	38-3240	38-3440
1.143	0.0450	2.819	0.111		38-3043	381-3043	381-3243	38-3243	38-3443
1.295	0.0510	2.972	0.117		38-3046	381-3046	381-3246	38-3246	38-3446
1.422	0.0560	3.099	0.112		38-3050	381-3050	381-3250	38-3250	38-3450
1.524	0.0600	3.200	0.126		38-3052	381-3052	381-3252	38-3252	38-3452
1.651	0.0650	3.327	0.131		38-3055	381-3055	381-3255	38-3255	38-3455
1.750	0.0690	3.430	0.135		38-3058	381-3058	381-3258	38-3258	38-3458
1.854	0.0730	3.531	0.139		38-3060	381-3060	381-3260	38-3260	38-3460
2.057	0.0810	3.734	0.147		38-3065	381-3065	381-3265	38-3265	38-3465
2.286	0.0900	3.962	0.156		38-3067	381-3067	381-3267	38-3267	38-3467
2.540	0.1000	4.216	0.166		38-3070	381-3070	381-3270	38-3270	38-3470
2.794	0.1100	4.470	0.176		38-3078	381-3078	381-3278	38-3278	38-3478

Silicone – by the meter

					By the meter  3 m (118.11 inch) Per package 25.00 €	By the meter  10 m (393.70 inch) Per package 50.00 €
Inside diameter:		Outside diameter:		Colour code:		
mm	inch	mm	inch		Part number	Part number
0.508	0.0200	2.337	0.092		38-3125	38-3525
0.635	0.0250	2.438	0.097		38-3130	38-3530
0.762	0.0300	2.464	0.096		38-3132	38-3532
0.889	0.0350	2.565	0.101		38-3135	38-3535
1.016	0.0400	2.692	0.106		38-3140	38-3540
1.143	0.0450	2.819	0.111		38-3143	38-3543
1.295	0.0510	2.972	0.117		38-3146	38-3546
1.422	0.0560	3.099	0.112		38-3150	38-3550
1.524	0.0600	3.200	0.126		38-3152	38-3552
1.651	0.0650	3.327	0.131		38-3155	38-3555
1.750	0.0690	3.430	0.135		38-3158	38-3558
1.854	0.0730	3.531	0.139		38-3160	38-3560
2.057	0.0810	3.734	0.147		38-3165	38-3565
2.286	0.0900	3.962	0.156		38-3167	38-3567
2.540	0.1000	4.216	0.166		38-3170	38-3570
2.794	0.1100	4.470	0.176		38-3178	38-3578

Silicone double pump tubing



Pairs of tubing connected with Y-pieces on both sides extended with a tube of 25cm					 2 bridges 152 mm 250 mm 406 mm 250 mm per package 6 pieces 167.00 €
Inside diameter:		Outside diameter:		Colour code:	
mm	inch	mm	inch		Part number
1.016	0.0400	2.692	0.106		38-3340
1.143	0.0450	2.819	0.111		38-3343
1.295	0.0510	2.972	0.117		38-3346
1.422	0.0560	3.099	0.112		38-3350
1.524	0.0600	3.200	0.126		38-3352
1.651	0.0650	3.327	0.131		38-3355
1.750	0.0690	3.430	0.135		38-3358
1.854	0.0730	3.531	0.139		38-3360
2.057	0.0810	3.734	0.147		38-3365
2.286	0.0900	3.962	0.156		38-3367
2.540	0.1000	4.216	0.166		38-3370
2.794	0.1100	4.470	0.176		38-3378




Standard-Tubing-Connector



Polypropylene

Part number	Connection	Configuration	For tubes with an inside diameter	Price w.o. VAT
37-0200	Olive	Y-connector	0.90 – 1.80 mm	3.80 €
37-0202	Olive	Y-connector	1.80 – 2.80 mm	3.80 €
37-0204	Olive	Y-connector	2.80 – 3.20 mm	3.80 €
37-0201	Olive	I-connector	0.90 – 1.80 mm	2.50 €
37-0203	Olive	I-connector	1.80 – 2.80 mm	2.50 €
37-0205	Olive	I-connector	2.80 – 3.20 mm	2.50 €
37-0208	Olive	I-reduction	0.90 – 1.80 to 1.80 – 2.80 mm	2.50 €
37-0210	Olive	I-reduction	1.80 – 2.80 to 2.80 – 3.20 mm	2.50 €
37-0211	Olive	I-reduction	0.90 – 1.80 to 2.80 – 3.20 mm	2.50 €

We are also happy to manufacture tubing already pre-fitted with tubing connectors.
Can't find the right tubing connector? We also offer Luer and other connectors. Please feel free to contact us.



Special tubing

Three simple steps to your custom tubing

1

Define your specification

Length	The dimensions of your tubing
Bridges	With or without bridges, any number
Bridge interval	The dimensions of your tubing
Material	PVC Standard, Solvent Flex, PU Longlife, fluoro rubber, Santoprene, silicone
Internal diameter	0.127 through 3.175 mm
Number of tubes per pack	Custom

Do you need a specific color and Shore hardness? Feel free to get in touch with us!

2

You will receive a personalized offer

We send out offers within 24 hours.

If the offer meets with your approval, you can order directly by phone or email.

And if you are not happy with our offer? Then let us know. We will work with you to find an ideal implementation and adjust the offer accordingly.

3

You receive the goods you ordered

When we receive your order, we will quickly send you an order confirmation detailing the delivery date.


Have you not yet found the right tubing for your application?

We would be delighted to configure your custom tubing based on your particular specifications. Please send us the following specifications:

- Raw material
- Internal diameter
- Tubing length
- Number of bridges and
- Distance between the bridges

... and we will send you a personalized offer.

Just go to the link.



<https://www.spetec.de/en/products/laboratory-equipment/pump-tubing-tubes>

We will be happy to answer any questions you have. Just phone us on **+49 8122 95909-0** or mail us at spetec@spetec.de!



Single-element and multi-element standards

ICP standards

Spetec offers a wide range of **NIST-certified** ICP standards from Inorganic Ventures.



Possible applications:

- ICP & ICP-MS
- Ion chromatography
- Atomic absorption spectrometry (AAS)
- Water quality control
- Wet chemistry



<https://www.spetec.de/en/downloads-for-laboratory-equipment>

OEM and standalone pumps for laboratory use

Peristaltic pumps

Perimax/OEM (ml)



<https://www.spetec.de/en/downloads-for-laboratory-equipment>

Syringe pumps

SYMAX/OEM (nl – ml)







[https://www.spetec.de/en/products/
laboratory-equipment/pump-tubing-tubes](https://www.spetec.de/en/products/laboratory-equipment/pump-tubing-tubes)

SPETEC® GmbH
Am Kletthamer Feld 15
85435 Erding
Germany

Tel. +49 8122 95909-0
Fax. +49 8122 95909-55

Email: spetec@spetec.de
www.spetec.de

