



PU Longlife Solvent Resistant Flex

Inside Diameter:		Outside Diameter:		Colour Code:	2 bridges				3 bridges				
		mm	inch		mm	inch	Part Number	12 pieces per package	Price	Part Number	12 pieces per package	Price	
0.381	0.0150	2.209	0.087	orange/green	38-2020	12 pieces per package	46.00 €	381-2020	12 pieces per package	51.00 €	38-2420	12 pieces per package	51.00 €
0.508	0.0200	2.337	0.092	orange/yellow	38-2025	12 pieces per package	46.00 €	381-2025	12 pieces per package	51.00 €	38-2425	12 pieces per package	51.00 €
0.635	0.0250	2.464	0.097	orange/white	38-2030	12 pieces per package	46.00 €	381-2030	12 pieces per package	51.00 €	38-2430	12 pieces per package	51.00 €
0.762	0.0300	2.438	0.096	black/black	38-2032	12 pieces per package	46.00 €	381-2032	12 pieces per package	51.00 €	38-2432	12 pieces per package	51.00 €
0.889	0.0350	2.565	0.101	orange/orange	38-2035	12 pieces per package	46.00 €	381-2035	12 pieces per package	51.00 €	38-2435	12 pieces per package	51.00 €
1.016	0.0400	2.692	0.106	white/white	38-2040	12 pieces per package	46.00 €	381-2040	12 pieces per package	51.00 €	38-2440	12 pieces per package	51.00 €
1.143	0.0450	2.819	0.111	red/red	38-2043	12 pieces per package	46.00 €	381-2043	12 pieces per package	51.00 €	38-2443	12 pieces per package	51.00 €
1.295	0.0510	2.972	0.117	grey/grey	38-2046	12 pieces per package	46.00 €	381-2046	12 pieces per package	51.00 €	38-2446	12 pieces per package	51.00 €
1.422	0.0560	3.099	0.112	yellow/yellow	38-2050	12 pieces per package	46.00 €	381-2050	12 pieces per package	51.00 €	38-2450	12 pieces per package	51.00 €
1.524	0.0600	3.200	0.126	yellow/blue	38-2052	12 pieces per package	46.00 €	381-2052	12 pieces per package	51.00 €	38-2452	12 pieces per package	51.00 €
1.651	0.0650	3.327	0.131	blue/blue	38-2055	12 pieces per package	46.00 €	381-2055	12 pieces per package	51.00 €	38-2455	12 pieces per package	51.00 €
1.750	0.0690	3.430	0.135	blue/green	38-2058	12 pieces per package	46.00 €	381-2058	12 pieces per package	51.00 €	38-2458	12 pieces per package	51.00 €
1.854	0.0730	3.531	0.139	green/green	38-2060	12 pieces per package	46.00 €	381-2060	12 pieces per package	51.00 €	38-2460	12 pieces per package	51.00 €
2.057	0.0810	3.734	0.147	purple/purple	38-2065	12 pieces per package	46.00 €	381-2065	12 pieces per package	51.00 €	38-2465	12 pieces per package	51.00 €
2.286	0.0900	3.962	0.156	purple/black	38-2067	12 pieces per package	46.00 €	381-2067	12 pieces per package	51.00 €	38-2467	12 pieces per package	51.00 €
2.540	0.1000	4.216	0.166	purple/orange	38-2070	12 pieces per package	46.00 €	381-2070	12 pieces per package	51.00 €	38-2470	12 pieces per package	51.00 €
2.794	0.1100	4.470	0.176	purple/white	38-2078	12 pieces per package	46.00 €	381-2078	12 pieces per package	51.00 €	38-2478	12 pieces per package	51.00 €
3.175	0.1250	4.851	0.191	black/white	38-2080	12 pieces per package	46.00 €	381-2080	12 pieces per package	51.00 €	38-2480	12 pieces per package	51.00 €

PU Longlife Solvent Resistant Flex – by the meter

Inside Diameter:		Outside Diameter:		Colour Code:	Part Number	by the meter  3 m (118.11 inch) per package 26.00 €	Part Number	by the meter  10 m (393.70 inch) per package 52.00 €
mm	inch	mm	inch					
0.381	0.0150	2.209	0.087	orange/green	38-2120	38-2520	38-2520	
0.508	0.0200	2.337	0.092	orange/yellow	38-2125	38-2525	38-2525	
0.635	0.0250	2.464	0.097	orange/white	38-2130	38-2530	38-2530	
0.762	0.0300	2.438	0.096	black/black	38-2132	38-2532	38-2532	
0.889	0.0350	2.565	0.101	orange/orange	38-2135	38-2535	38-2535	
1.016	0.0400	2.692	0.106	white/white	38-2140	38-2540	38-2540	
1.143	0.0450	2.819	0.111	red/red	38-2143	38-2543	38-2543	
1.295	0.0510	2.972	0.117	grey/grey	38-2146	38-2546	38-2546	
1.422	0.0560	3.099	0.112	yellow/yellow	38-2150	38-2550	38-2550	
1.524	0.0600	3.200	0.126	yellow/blue	38-2152	38-2552	38-2552	
1.651	0.0650	3.327	0.131	blue/blue	38-2155	38-2555	38-2555	
1.750	0.0690	3.430	0.135	blue/green	38-2158	38-2558	38-2558	
1.854	0.0730	3.531	0.139	green/green	38-2160	38-2560	38-2560	
2.057	0.0810	3.734	0.147	purple/purple	38-2165	38-2565	38-2565	
2.286	0.0900	3.962	0.156	purple/black	38-2167	38-2567	38-2567	
2.540	0.1000	4.216	0.166	purple/orange	38-2170	38-2570	38-2570	
2.794	0.1100	4.470	0.176	purple/white	38-2178	38-2578	38-2578	
3.175	0.1250	4.851	0.191	black/white	38-2180	38-2580	38-2580	