


Solvent Flex – Lösemittelbeständig

Innendurchmesser:	Aussendurchmesser:		Farbcodierung:	2 Reiter				3 Reiter				
	mm	(Zoll)		mm	(Zoll)	je Packung	Bestellnummer	je Packung	Bestellnummer	je Packung	Bestellnummer	je Packung
0,381	0,0150	2,209	0,087	orange/grün	38-4020	381-4020	381-4220	38-4220	38-4420	38-4220	38-4420	38-4420
0,508	0,0200	2,337	0,092	orange/gelb	38-4025	381-4025	381-4225	38-4225	38-4425	38-4225	38-4425	38-4425
0,635	0,0250	2,464	0,097	orange/weiß	38-4030	381-4030	381-4230	38-4230	38-4430	38-4230	38-4430	38-4430
0,762	0,0300	2,438	0,096	schwarz/schwarz	38-4032	381-4032	381-4232	38-4232	38-4432	38-4232	38-4432	38-4432
0,889	0,0350	2,565	0,101	orange/orange	38-4035	381-4035	381-4235	38-4235	38-4435	38-4235	38-4435	38-4435
1,016	0,0400	2,692	0,106	weiß/weiß	38-4040	381-4040	381-4240	38-4240	38-4440	38-4240	38-4440	38-4440
1,143	0,0450	2,819	0,111	rot/rot	38-4043	381-4043	381-4243	38-4243	38-4443	38-4243	38-4443	38-4443
1,295	0,0510	2,972	0,117	grau/grau	38-4046	381-4046	381-4246	38-4246	38-4446	38-4246	38-4446	38-4446
1,422	0,0560	3,099	0,112	gelb/gelb	38-4050	381-4050	381-4250	38-4250	38-4450	38-4250	38-4450	38-4450
1,524	0,0600	3,200	0,126	gelb/blau	38-4052	381-4052	381-4252	38-4252	38-4452	38-4252	38-4452	38-4452
1,651	0,0650	3,327	0,131	blau/blau	38-4055	381-4055	381-4255	38-4255	38-4455	38-4255	38-4455	38-4455
1,750	0,0690	3,430	0,135	blau/grün	38-4058	381-4058	381-4258	38-4258	38-4458	38-4258	38-4458	38-4458
1,854	0,0730	3,531	0,139	grün/grün	38-4060	381-4060	381-4260	38-4260	38-4460	38-4260	38-4460	38-4460
2,057	0,0810	3,734	0,147	violett/violett	38-4065	381-4065	381-4265	38-4265	38-4465	38-4265	38-4465	38-4465
2,286	0,0900	3,962	0,156	violett/schwarz	38-4067	381-4067	381-4267	38-4267	38-4467	38-4267	38-4467	38-4467
2,540	0,1000	4,216	0,166	violett/orange	38-4070	381-4070	381-4270	38-4270	38-4470	38-4270	38-4470	38-4470
2,794	0,1100	4,470	0,176	violett/weiß	38-4078	381-4078	381-4278	38-4278	38-4478	38-4278	38-4478	38-4478
3,175	0,1250	4,851	0,191	schwarz/weiß	38-4080	381-4080	381-4280	38-4280	38-4480	38-4280	38-4480	38-4480

Solvent Flex – Lösemittelbeständig – Meterware

Innendurchmesser:		Aussendurchmesser:		Farbcodierung:	Bestellnummer	Meterware je Packung 3 m Stück 21,00 € zzgl. MwSt.
mm	(Zoll)	mm	(Zoll)			
0,381	0,0150	2,209	0,087	orange/grün	38-4120	 je Packung 10 m Stück 42,00 € zzgl. MwSt.
0,508	0,0200	2,337	0,092	orange/gelb	38-4125	
0,635	0,0250	2,464	0,097	orange/weiß	38-4130	
0,762	0,0300	2,438	0,096	schwarz/schwarz	38-4132	
0,889	0,0350	2,565	0,101	orange/orange	38-4135	
1,016	0,0400	2,692	0,106	weiß/weiß	38-4140	
1,143	0,0450	2,819	0,111	rot/rot	38-4143	
1,295	0,0510	2,972	0,117	grau/grau	38-4146	
1,422	0,0560	3,099	0,112	gelb/gelb	38-4150	
1,524	0,0600	3,200	0,126	gelb/blau	38-4152	
1,651	0,0650	3,327	0,131	blau/blau	38-4155	
1,750	0,0690	3,430	0,135	blau/grün	38-4158	
1,854	0,0730	3,531	0,139	grün/grün	38-4160	
2,057	0,0810	3,734	0,147	violett/violett	38-4165	
2,286	0,0900	3,962	0,156	violett/schwarz	38-4167	
2,540	0,1000	4,216	0,166	violett/orange	38-4170	
2,794	0,1100	4,470	0,176	violett/weiß	38-4178	
3,175	0,1250	4,851	0,191	schwarz/weiß	38-4180	