

























































Solvent Flex

Inside Diameter:		Outside Diameter:		Colour Code:	2 bridges				3 bridges							
		mm	inch		mm	inch	Part Number	12 pieces per package	Price	Part Number	12 pieces per package	Price				
0.381	0.0150	2.209	0.087	orange/green	38-4020	12 pieces per package	37.00 €		381-4020	12 pieces per package	41.00 €		381-4220	12 pieces per package	41.00 €	
0.508	0.0200	2.337	0.092	orange/yellow	38-4025	12 pieces per package	37.00 €		381-4025	12 pieces per package	41.00 €		381-4225	12 pieces per package	41.00 €	
0.635	0.0250	2.464	0.097	orange/white	38-4030	12 pieces per package	37.00 €		381-4030	12 pieces per package	41.00 €		381-4230	12 pieces per package	41.00 €	
0.762	0.0300	2.438	0.096	black/black	38-4032	12 pieces per package	37.00 €		381-4032	12 pieces per package	41.00 €		381-4232	12 pieces per package	41.00 €	
0.889	0.0350	2.565	0.101	orange/orange	38-4035	12 pieces per package	37.00 €		381-4035	12 pieces per package	41.00 €		381-4235	12 pieces per package	41.00 €	
1.016	0.0400	2.692	0.106	white/white	38-4040	12 pieces per package	37.00 €		381-4040	12 pieces per package	41.00 €		381-4240	12 pieces per package	41.00 €	
1.143	0.0450	2.819	0.111	red/red	38-4043	12 pieces per package	37.00 €		381-4043	12 pieces per package	41.00 €		381-4243	12 pieces per package	41.00 €	
1.295	0.0510	2.972	0.117	grey/grey	38-4046	12 pieces per package	37.00 €		381-4046	12 pieces per package	41.00 €		381-4246	12 pieces per package	41.00 €	
1.422	0.0560	3.099	0.112	yellow/yellow	38-4050	12 pieces per package	37.00 €		381-4050	12 pieces per package	41.00 €		381-4250	12 pieces per package	41.00 €	
1.524	0.0600	3.200	0.126	yellow/blue	38-4052	12 pieces per package	37.00 €		381-4052	12 pieces per package	41.00 €		381-4252	12 pieces per package	41.00 €	
1.651	0.0650	3.327	0.131	blue/blue	38-4055	12 pieces per package	37.00 €		381-4055	12 pieces per package	41.00 €		381-4255	12 pieces per package	41.00 €	
1.750	0.0690	3.430	0.135	blue/green	38-4058	12 pieces per package	37.00 €		381-4058	12 pieces per package	41.00 €		381-4258	12 pieces per package	41.00 €	
1.854	0.0730	3.531	0.139	green/green	38-4060	12 pieces per package	37.00 €		381-4060	12 pieces per package	41.00 €		381-4260	12 pieces per package	41.00 €	
2.057	0.0810	3.734	0.147	purple/purple	38-4065	12 pieces per package	37.00 €		381-4065	12 pieces per package	41.00 €		381-4265	12 pieces per package	41.00 €	
2.286	0.0900	3.962	0.156	purple/black	38-4067	12 pieces per package	37.00 €		381-4067	12 pieces per package	41.00 €		381-4267	12 pieces per package	41.00 €	
2.540	0.1000	4.216	0.166	purple/orange	38-4070	12 pieces per package	37.00 €		381-4070	12 pieces per package	41.00 €		381-4270	12 pieces per package	41.00 €	
2.794	0.1100	4.470	0.176	purple/white	38-4078	12 pieces per package	37.00 €		381-4078	12 pieces per package	41.00 €		381-4278	12 pieces per package	41.00 €	
3.175	0.1250	4.851	0.191	black/white	38-4080	12 pieces per package	37.00 €		381-4080	12 pieces per package	41.00 €		381-4280	12 pieces per package	41.00 €	

Solvent Flex – by the meter

Inside Diameter:		Outside Diameter:		Colour Code:	Part Number	by the meter  3 m (118.11 inch) per package 21.00 €	Part Number	by the meter  10 m (393.70 inch) per package 42.00 €
mm	inch	mm	inch					
0.381	0.0150	2.209	0.087	orange/green	38-4120	38-4520	38-4520	
0.508	0.0200	2.337	0.092	orange/yellow	38-4125	38-4525	38-4525	
0.635	0.0250	2.464	0.097	orange/white	38-4130	38-4530	38-4530	
0.762	0.0300	2.438	0.096	black/black	38-4132	38-4532	38-4532	
0.889	0.0350	2.565	0.101	orange/orange	38-4135	38-4535	38-4535	
1.016	0.0400	2.692	0.106	white/white	38-4140	38-4540	38-4540	
1.143	0.0450	2.819	0.111	red/red	38-4143	38-4543	38-4543	
1.295	0.0510	2.972	0.117	grey/grey	38-4146	38-4546	38-4546	
1.422	0.0560	3.099	0.112	yellow/yellow	38-4150	38-4550	38-4550	
1.524	0.0600	3.200	0.126	yellow/blue	38-4152	38-4552	38-4552	
1.651	0.0650	3.327	0.131	blue/blue	38-4155	38-4555	38-4555	
1.750	0.0690	3.430	0.135	blue/green	38-4158	38-4558	38-4558	
1.854	0.0730	3.531	0.139	green/green	38-4160	38-4560	38-4560	
2.057	0.0810	3.734	0.147	purple/purple	38-4165	38-4565	38-4565	
2.286	0.0900	3.962	0.156	purple/black	38-4167	38-4567	38-4567	
2.540	0.1000	4.216	0.166	purple/orange	38-4170	38-4570	38-4570	
2.794	0.1100	4.470	0.176	purple/white	38-4178	38-4578	38-4578	
3.175	0.1250	4.851	0.191	black/white	38-4180	38-4580	38-4580	